

## Fixing and Staining Worms on a Slide

1. Use 1%, 1.5%, and 2% PFA with a 5-minute fix. Determine empirically which one works best for your antibody.
  - a. Place worms in 5 $\mu$ l 1X sperm salts on a positively-charged glass slide
  - b. Using a needle (Precision Glide<sup>®</sup> Needle, Becton Dickinson & Co., 22 G 1  $\frac{1}{2}$ ), make an incision between the pharynx and gonad (closer to the gonad than the pharynx) or between the tip of the tail and the gonad of each worm. To release the embryos, cut at the vulva.
  - c. Add 5  $\square$ l of 2%, 3%, or 4% PFA (Electron Microscopy Sciences, cat. #15710) in 1X sperm salts, mix a little.

To prepare 2% PFA:	125 $\mu$ l 16% PFA 200 $\mu$ l 5X sperm salts 675 $\mu$ l water
To prepare 3% PFA:	187.5 $\mu$ l 16% PFA 200 $\mu$ l 5X sperm salts 612.5 $\mu$ l water
To prepare 4% PFA:	200 $\mu$ l 16% PFA 200 $\mu$ l 5X sperm salts 550 $\mu$ l water
  - d. Place the slide in a humid chamber for 5 minutes without coverslip. Place an 18-mm x 18-mm coverslip on top of the worms. Freeze the slide on a piece of dry ice. Let the slide sit on the dry ice for at least 10 minutes.
2. Remove the slide from the dry ice and quickly pop off the coverslip with a quick downward stroke of a single razor.
3. Optional: Immediately place the slide in 95% ethanol for 1 minutes. Be careful or the worms will fall off. (Ethanol coagulates proteins – can do without this step.)
4. Incubate the slide in a vessel containing PBST for at least 10 minutes (Longer washes can decrease background but reverse fix). Repeat this step twice.
5. Wick off most of the PBST, but do not allow the worms to dry. Add 30  $\mu$ l primary antibody diluted in PBST. Carefully cover the worms with a 20 x 20-mm piece of parafilm. Place the slide in a humid chamber. Incubate overnight at room temperature (or for a few hours at 37°C).
6. Incubate the slide in a vessel containing PBST for at least ten minutes. Carefully remove the piece of parafilm.

Adapted from Annette Chan  
January 15, 1999

7. Incubate the slide in a vessel containing fresh PBST for at least 10 minutes. Repeat this step once.
8. Wick off most of the PBST, but do not allow the worms to dry. Add 25  $\mu$ l secondary antibody diluted in PBST. Carefully cover the worms with a 20 x 20-mm piece of parafilm. Place the slide in a humid chamber. Incubate for 4 hours at room temperature (30 min to 1 hr at 37°C). Keep the slides in the dark whenever possible.
9. Incubate the slide in a vessel containing PBST for at least 10 minutes. Carefully remove the piece of parafilm.
10. Incubate the slide in a vessel containing fresh PBST for at least ten minutes. Repeat this step with fresh PBST to which you add 1  $\mu$ l of 1 mg/mL DAPI.
11. Wick off most of the PBST, but do not allow the worms to dry.
12. Mount the slide with 20  $\mu$ l Vectashield or other mounting medium and a 22 x 40-mm cover glass. Seal the slide with fingernail polish. Alternatively, use a 22 X 50 mm coverglass and do not seal the slides.
13. Store the slides at -20°C. Keep the slides dark whenever possible.

<b>1X Sperm Salts:</b>	50 mM PIPES, pH7.0	<b>5X:</b>	250 mM PIPES pH 7.0
	25 mM KCl		125 mM KCl
	1 mM MgSO <sub>4</sub>		5 mM MgSO <sub>4</sub>
	45 mM NaCl		225 mM NaCl
	2 mM CaCl <sub>2</sub>		10 mM CaCl <sub>2</sub>

**PBST:** 100 ml 10X PBS  
5 ml Triton X-100  
2 ml 0.5 M EDTA, pH 8  
ddH<sub>2</sub>O to 1L

**DABCO (Mounting medium):**

5 g DABCO  
5 ml PBS  
45 ml glycerin