

Cell Cycle Analysis of *Drosophila* Kc167 cells using Alexa 488 Click-IT EdU and FxCycle Violet on the Attune Cytometer

Dec. 14th 2011

Laura Buttitta
University of Michigan, MCDB Dept.

Drosophila Kc167 cells (<https://dgrc.cgb.indiana.edu/product/View?product=1>) are cultured in complete medium at 25°C. (Schneider's Medium Gibco #21720-024, 10% heat-inactivated FBS HyClone SH30070.03, 1X Pen-Strep Gibco #15070-063).

Cells were plated at a density of 1×10^6 cells/mL in a 6-well plate in complete medium containing either DMSO (control) or 20 μ M SP600125 (Sigma S5567) to arrest cells, for 48 hours. (Note that in contrast to mammalian cells, proliferating *Drosophila* cells in complete medium are primarily in G2 phase¹).

Cells were exposed to 10 μ M EdU at 25°C in complete medium for 30 min. Cells were re-suspended by gentle pipetting and transferred to an eppendorf (*Drosophila* Kc167 cells are semi-adherent).

Staining Protocol:

- Pellet cells by centrifugation at 1,000K RPM room temp. for 5 min., and discard media.
- Wash cells with 1mL 1XPBX+1%BSA
- Pellet cells at 1000K RPM room temp. for 5 min. and discard the wash.
- Resuspend cells in 50 μ L of Click-IT fixative, incubate at room temp. for 15 min.
- Pellet cells at 1000K RPM room temp. for 5 min. and discard the fix.
- Resuspend pellet in 1mL 1XPBX+1%BSA to wash.
- Pellet cells at 1000K RPM room temp. for 5 min. and discard the wash.
- Resuspend pellet in a second wash of 1mL 1XPBX+1%BSA.
- Pellet cells again at 1000K RPM room temp. for 5 min. and discard the 2nd wash.
- Resuspend pellet in 30 μ L of 1X Saponin Permeabilization buffer.
- Immediately prepare Click-IT EdU cocktail: (use 150 μ L cocktail per sample, for 500 μ L cocktail prepare in order):
 - 438 μ L 1XPBS
 - 10 μ L CuSO₄
 - 5 μ L 488-Click-IT Azide (this is 2X the standard protocol, our lab stock is diluted 2X)
 - 50 μ L 1X Buffer additive (this must be freshly prepared from 10X stock by dilution with dH₂O and used within 15 min. of dilution)
- Mix well, add 150 μ L cocktail per sample. –Must be used within 15 min. of preparation.

- Incubate for 30 min. room temp. protected from light.
 - Pellet cells at 1000K RPM room temp. for 5 min. and discard supernatant.
 - Wash pellet in 100 μ L of 1X Saponin Permeabilization buffer.
 - Pellet cells again at 1000K RPM room temp. for 5 min. and discard wash.
 - Resuspend pellet in 100 μ L of 1X Saponin Permeabilization buffer.
 - Dilute FxCycle Violet in 1:1,000 in 1XPBS. Add 1mL FxCycle Violet in 1XPBS to each sample in 100 μ L of 1X Saponin Permeabilization buffer.
 - Incubate at room temp. for 30 min.
-
- Run cells on the Attune cytometer using standard sensitivity at a flow rate of 100uL/min. (Cells can be run on cytometer directly from the eppendorf.)